

Autoclaves

72-200 Liter Chamber



Autoclaves for laboratories and hospitals

Compact floor device with a rectangular or round sterilization chamber and automatic locking door for easy handling.

HST-Series	3x3x6	50/75	50/100
Effective space (mm)	320 x 350 x 630	Ø 500 x 750	Ø 500 x 1000
Chamber volume (liter)	72	150	200
Outer dimensions (mm)	600 x 1450 x 850	800 x 1800 x 955	800 x 1800 x 1205

Excerpt from the extensive special options:

- Pass-through version with gas-tight separation (Bioseal)
- Exhaust filter with condensate inactivation, for the sterilization of infectious goods of risk-Cl. S2 and S3
- Vacuum pump set with water saving device for single-stage or fractionated pre-vacuum and final drying under vacuum
- Transport and loading trolley for ergonomic loading and unloading.

