

Sterilizer HST 4x6x6

Floor device

Compact floor device with a rectangular sterilization chamber and full automatic vertical sliding door for easy handling. Small overall dimensions with a high loading capacity take care for an efficient work.

Possible accessory:

- Rapid water recooling together with the supporting pressure regulation, this ensures the shortest re-cooling times of liquid media via a sterile air filter.
- Exhaust air filter with condensate inactivation for the sterilization of infectious materials of risk class S2 and S3.
Filter retention rate: > 99.5% 0.02µm
- Vacuum pump set, to safely evacuate air from the sterilized items using a single-stage or fractionated pre-vacuum. Enables the sterilized material to be removed when dry at the end of the process.
- Transport and loading trolleys for ergonomic loading

Technical data

Effective space mm	Usable chamber volume	Load height	Overall dimension mm
(w)450 x (h)700 x (d)650	205 Liter	780 mm	(w)800 x (h)1950 x (d)1000

