

# Laboratory Tabletop Autoclave

24-72 Litres Volume

The laboratory tabletop autoclaves distinguish particularly through their big capacity compared with their compact external size.

The LTA 3x3x6 has due to the squared chamber a capacity up to 4 pieces of 5-litres lab bottles!

LTA-Serie	275	2x3x4	3x3x6
Usable space (mm)	Ø 294 x 350	280 x 320 x 420	320 x 350 x 630
Chamber (litres)	24	42	73
External dimension (mm)	500 x 650 x 670		700 x 800 x 750



LTA 3x3x6  
Squared Chamber

Completely made of stainless steel.

Extract of the extensive options available:

- Magnetic mixer system for still shorter process times
- Batch documentation via interface or USB-stick.
- Water quick re-cooling „FAST-COOLER“, heat exchanger, fan within.



LTA 275  
Round Chamber