

Laboratory vertical autoclave



Easy to operate thanks to a touchscreen with a high program variability

Short cycles with a low media consumption

Robust and long-lasting

Extract of the extensive options:

- Built-in steam generator Independent of the sterilisation chamber
- Water quick re-cooling „FASTCOOLER“ Heat exchanger, fan within.
- Vacuum pump assembly for single-step or fractionated prevacuum and drying under vacuum
- Exhaust air filter with condensate inactivation, for the sterilisation of infectious goods of the risk class S2

LVSA-Serie	40/60	40/75	50/60	50/70
Usable space in mm	(Ø)400 x (H)620	(Ø)400 x (H)770	(Ø)500x (H)680	(Ø)500x (H)780
Chamber volume in Liter	80	96	133	153
External dimension in mm	(B) 550 x (T) 750		(B) 750 x (T) 900	

